

Vision Cam 1000

Affordable High-Resolution Live Imaging Camera.

Description

The Vision Cam 1000 is a competitively priced high-resolution live imaging microscopy camera. With a live image of 1280x1024 pixels, this image can be projected in any lecture or classroom environment without fear of pixelation. Enjoy crisp and clear live images.

Images Plus is an imaging suite included for both PC and Macintosh OSX platforms providing measurement, analysis and sharing tools for teachers, students and researchers alike.

With all necessary accessories included, the Vision Cam 1000 fits virtually any Microscope's eyepiece therefore allowing you to turn any microscope into a Digital Microscope.



Specifications

Camera & Accessories	
Chip Configuration	1/3" (optical) Live 1.3 Mega Pixel Imaging Chip in portable casing
Data Output	USB2.0
Power Supply	Through USB Cable from Computer
Included Software	Images Plus 2.0ML for Windows & Images Plus 2.0 for Macintosh OSX
Included Accessories	12mm Focusable Lens 28, 30, 34 & 35mm Eyepiece Couplers Macro Tube Calibration Slide for accurate measurements
Min. Specifications	Windows: XP or Higher; 256MB RAM; P4; Built-in USB2.0 Macintosh: OSX or higher; 256MB RAM; G4; Built-in USB2.0



Vision Engineering manufacture a comprehensive range of ergonomic stereo microscopes as well as a complete line of optical and video non-contact measuring systems.

For more information...

Vision Engineering has a network of offices and technical distributors around the world. For more information, please contact your Vision Engineering branch, local authorised distributor, or visit our website.

Vision Engineering Ltd.
(Manufacturing)
Send Road, Send,
Woking, Surrey, GU23 7ER, England
Tel: +44 (0) 1483 248300
Email: generalinfo@visioneng.com

Vision Engineering Ltd.
(Commercial)
Monument House, Monument Way West,
Woking, Surrey, GU21 5EN, England
Tel: +44 (0) 1483 248300
Email: generalinfo@visioneng.com

Vision Engineering Inc.
(Manufacturing & Commercial)
570 Danbury Road, New Milford,
CT 06776 USA
Tel: +1 (860) 355 3776
Email: info@visioneng.com

Vision Engineering Inc.
(West Coast Commercial)
745 West Taft Avenue, Orange,
CA 92865 USA
Tel: +1 (714) 974 6966
Email: info@visioneng.com

Vision Engineering Ltd.
(Central Europe)
Anton-Pendele-Str. 3,
82275 Emmering, Deutschland
Tel: +49 (0) 8141 40167-0
Email: info@visioneng.de

Vision Engineering Ltd.
(France)
ZAC de la Tremblaie, Av. de la Tremblaie
91220 Le Plessis Paté, France
Tel: +33 (0) 160 76 60 00
Email: info@visioneng.fr

Vision Engineering Ltd.
(Italia)
Via Cesare Cantù, 9
20092 Cinisello Balsamo MI, Italia
Tel: +39 02 6129 3518
Email: info@visioneng.it

Nippon Vision Engineering
(Japan)
272-2 Saedo-cho, Tsuduki-ku,
Yokohama-shi, 224-0054, Japan
Tel: +81 (0) 45 935 1117
Email: info@visioneng.jp

Vision Engineering Ltd
(China)
11J, International Ocean Building,
720 Pudong Avenue, Shanghai,
200120, P.R. China
Tel: +86 (0) 21 5036 7556
Email: info@visioneng.com.cn

Vision Engineering
(S.E. Asia)
Tel: +603 80700908
Email: info@visioneng.asia

Vision Engineering
(India)
Email: info@visioneng.co.in

Distributor



Visit our multi-lingual website:

www.visioneng.com